



Evaluation of the impact of bad storage (temperature) on the chicken meat quality by Raman Spectroscopy

A. ASSAF, E. GRANGÉ, S. JOUANNEAU, J. LAHMAR, M.J. DURAND and G. THOUAND

University of Nantes, UMR CNRS GEPEA 6144, TEAM-LRSY,
18 Boulevard Gaston Defferre, 85035 La Roche sur Yon, France

REAMIT Networking Symposium

9th January 2020
Crowne Plaza Nottingham, UK




Nantes University: A leading of higher education & research center in Western France

38 000 students

64 Research Units

3 754 staff members (2 118 assistant, associate and full professors and researchers + 1 636 administrative and technical staff)



GEPEA: process engineering for environment and food


220 members (associate and full professors, researchers + administrative and technical staff)

5 research topics


Bioprocess applied
on microalgae




Matrices & Food :
Processes /
Properties /
Structure - Sensory



Water Treatment-
Air- Metrology



Matrices &
Optimization -
System - Energy



Energy/Residuals
and Emission
Treatment



GEPEA: process engineering for environment and food

220 members (associate and full professors, researchers + administrative and technical staff)

5 research topics


Bioprocess applied
on microalgae




Matrices & Food :
Processes /
Properties /
Structure - Sensory



Water Treatment-
Air- Metrology



Matrices &
Optimization -
System - Energy




Energy/Residuals
and Emission
Treatment



GEPEA: TEAM-LRSY


Optical methods
Raman,
Fluorescence, ...



Different field of applications


GEPEA: TEAM-LRSY

Optical methods
Raman,
Fluorescence, ...




- Pollution fingerprints, toxicity, biodegradation**
- Microalgae process monitoring**
- Microorganism's detection**
- Food characterization**
- + other applications**

Raman scattering ($1/10^9$ photons)




Raman spectra




λ (Laser)


+ other radiative emissions or not
(Fluorescence, phosphorescence,)



Raman scattering ($1/10^9$ photons)




+ other radiative emissions or not
(Fluorescence, phosphorescence,)



North-West Europe


REA
European Region





- Non-invasive and non-destructive analysis:
no sample preparation required
- A **small quantity** of sample is sufficient
- **Very fast** analysis
- One shot = **simultaneous verification**
The presence/absence of pathogens
Lipid/protein/carbohydrate,...



NON DESTRUCTIVE



3D ANALYSIS




Optimization of optical parameters

- Laser: λ , P, t, R, ...
- Biology: high variability of living samples
- Statistics: nb of spectra, data, model reliability


All of these parameter should be examined to provide reliable data

Selection of the best wavelength



Raman spectra done on chicken matrices with 2 different wavelength (A) 532 nm & (B) 785 nm.


Variability between batches



Correlation between 2 batches of chicken


North-West Europe



REAMIT



North-West Europe


Raman spectra of chicken samples (left side) and 2D plots of spectra showing the Raman shift and the intensity differences between the studied samples (right side)





Kruskal-Wallis ANOVA Table

Source	SS	df	MS	Chi-sq	Prob>Chi-sq
Groups	221917.7	4	55479.4	117.57	1.7624e-24
Error	59319.8	145	409.1		
Total	281237.5	149			




J0



J4

North-West Europe

Classification of the Raman spectra done on chicken samples according to their storage conditions




North-West Europe

REAMIT

European Regional Development Fund

UoN will perform the following activities in REAMIT

- (i) Raman Spectroscopy technology for monitoring food quality in the supply chain
- (ii) Send the sensor data to Big data platforms for further analytics
- (iii) Conduct technology testing and piloting
- (iv) Help in project evaluation, risk assessment, monitoring and dissemination.



Raman spectroscopy is a powerful technique to monitor the quality of food matrices

Challenges for REAMIT project



Thank you for your attention

